



UFM1 Antibody, Biotin conjugated

Product Code	CSB-PA025560LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61960
Immunogen	Recombinant Human Ubiquitin-fold modifier 1 protein (11-63AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin-fold modifier 1, UFM1, C13orf20
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	UFM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.