



USP6 Antibody, Biotin conjugated

Product Code	CSB-PA025747LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35125
Immunogen	Recombinant Human Ubiquitin carboxyl-terminal hydrolase 6 protein (1122-1359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquitin carboxyl-terminal hydrolase 6 (EC 3.4.19.12) (Deubiquitinating enzyme 6) (Proto-oncogene TRE-2) (Ubiquitin thioesterase 6) (Ubiquitin-specific-processing protease 6), USP6, HRP1 TRE2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	USP6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.