



XPO1 Antibody, Biotin conjugated

Product Code	CSB-PA026221LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14980
Immunogen	Recombinant Human Exportin-1 protein (625-845AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Exportin-1 (Exp1) (Chromosome region maintenance 1 protein homolog), XPO1, CRM1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	XPO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.