



GPX4 Antibody, Biotin conjugated

Product Code	CSB-PA03049D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36969
Immunogen	Recombinant Human Phospholipid hydroperoxide glutathione peroxidase, mitochondrial protein (74-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipid hydroperoxide glutathione peroxidase (PHGPx) (EC 1.11.1.12) (Glutathione peroxidase 4) (GPx-4) (GSHPx-4), GPX4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPX4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.