



HSPA8 Antibody, Biotin conjugated

Product Code	CSB-PA03249D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11142
Immunogen	Recombinant Human Heat shock cognate 71 kDa protein (418-597AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8) (Lipopolysaccharide-associated protein 1) (LAP-1) (LPS-associated protein 1), HSPA8, HSC70 HSP73 HSPA10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HSPA8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.