



DUOXA1 Antibody, FITC conjugated

Product Code	CSB-PA04949C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q1HG43
Immunogen	Recombinant Human Dual oxidase maturation factor 1 protein (73-183AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual oxidase maturation factor 1 (Dual oxidase activator 1) (Numb-interacting protein), DUOXA1, NIP NUMBIP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DUOXA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.