



RPL13A Antibody, Biotin conjugated

Product Code	CSB-PA05289D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40429
Immunogen	Recombinant Human 60S ribosomal protein L13a protein (63-190AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L13a (23 kDa highly basic protein) (Large ribosomal subunit protein uL13), RPL13A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL13A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.