



GSK3A Antibody, Biotin conjugated

Product Code	CSB-PA05369D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49840
Immunogen	Recombinant Human Glycogen synthase kinase-3 alpha protein (36-118AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycogen synthase kinase-3 alpha (GSK-3 alpha) (EC 2.7.11.26) (Serine/threonine-protein kinase GSK3A) (EC 2.7.11.1), GSK3A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GSK3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.