



CXCL2 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA08989C0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P19875 |
| Immunogen | Recombinant Human C-X-C motif chemokine 2 protein (35-107AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | C-X-C motif chemokine 2 (Growth-regulated protein beta) (Gro-beta) (Macrophage inflammatory protein 2-alpha) (MIP2-alpha) [Cleaved into: GRO-beta(5-73) (GRO-beta-T) (Hematopoietic synergistic factor) (HSF) (SB-251353)], CXCL2, GRO2 GROB MIP2A SCYB2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CXCL2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |