



CXCL2 Antibody, Biotin conjugated

Product Code	CSB-PA08989D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19875
Immunogen	Recombinant Human C-X-C motif chemokine 2 protein (35-107AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-X-C motif chemokine 2 (Growth-regulated protein beta) (Gro-beta) (Macrophage inflammatory protein 2-alpha) (MIP2-alpha) [Cleaved into: GRO-beta(5-73) (GRO-beta-T) (Hematopoietic synergistic factor) (HSF) (SB-251353)], CXCL2, GRO2 GROB MIP2A SCYB2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CXCL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.