



CD247 Antibody, Biotin conjugated

Product Code	CSB-PA10499D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20963
Immunogen	Recombinant Human T-cell surface glycoprotein CD3 zeta chain protein (88-153AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell surface glycoprotein CD3 zeta chain (T-cell receptor T3 zeta chain) (CD antigen CD247), CD247, CD3Z T3Z TCRZ
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD247
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.