




IL17RA Antibody

Product Code	CSB-PA12929A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F46
Immunogen	The antiserum was produced against synthesized peptide derived from human IL17RA. AA range:141-190
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL17RA

Image



Western Blot

Positive WB detected in: HL60 whole cell lysate(30µg), THP-1 whole cell lysate(30µg), PC-3 whole cell lysate(30µg)

All lanes: IL17RA antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution

Predicted band size: 97,93 kDa

Observed band size: 97 kDa

Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.