



IL17RA Antibody, Biotin conjugated

Product Code	CSB-PA12929D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96F46
Immunogen	Recombinant Human Interleukin-17 receptor A protein (552-743AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-17 receptor A (IL-17 receptor A) (IL-17RA) (CDw217) (CD antigen CD217), IL17RA, IL17R
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL17RA