



FGF23 Antibody, HRP conjugated

Product Code	CSB-PA13919B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZV9
Immunogen	Recombinant Human Fibroblast growth factor 23 protein (123-221AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 23 (FGF-23) (Phosphatonin) (Tumor-derived hypophosphatemia-inducing factor) [Cleaved into: Fibroblast growth factor 23 N-terminal peptide; Fibroblast growth factor 23 C-terminal peptide], FGF23, HYPF
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGF23
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.