



PRKCG Antibody, Biotin conjugated

Product Code	CSB-PA15869D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05129
Immunogen	Recombinant Human Protein kinase C gamma type protein (1-314AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein kinase C gamma type (PKC-gamma) (EC 2.7.11.13), PRKCG, PKCG
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PRKCG
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.