



DST Antibody, FITC conjugated

Product Code	CSB-PA17019C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03001
Immunogen	Recombinant Human Dystonin protein (1-195AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dystonin (230 kDa bullous pemphigoid antigen) (230/240 kDa bullous pemphigoid antigen) (Bullous pemphigoid antigen 1) (BPA) (Bullous pemphigoid antigen) (Dystonia musculorum protein) (Hemidesmosomal plaque protein), DST, BP230 BP240 BPAG1 DMH DT KIAA0728
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DST
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.