



ATP1A3 Antibody

Product Code	CSB-PA18139A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13637
Immunogen	Recombinant Human Sodium/potassium-transporting ATPase subunit alpha-3 protein (143-278AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium/potassium-transporting ATPase subunit alpha-3 (Na(+)/K(+) ATPase alpha-3 subunit) (EC 3.6.3.9) (Na(+)/K(+) ATPase alpha(III) subunit) (Sodium pump subunit alpha-3), ATP1A3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ATP1A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.