



Allergen Fel d 4 Antibody, Biotin conjugated

Product Code	CSB-PA19827D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VFH6
Immunogen	Recombinant Cat Allergen Fel d 4 protein (16-186AA)
Raised In	Rabbit
Species Reactivity	Cat
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Allergen Fel d 4 (allergen Fel d 4)
Immunogen Species	Felis catus (Cat) (Felis silvestris catus)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.