



Major allergen Bla g 1.0101 Antibody, Biotin conjugated

Product Code	CSB-PA19837D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UAM5
Immunogen	Recombinant Blattella germanica Major allergen Bla g 1.0101 protein (1-412AA)
Raised In	Rabbit
Species Reactivity	Blattella germanica
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major allergen Bla g 1.0101
Immunogen Species	Blattella germanica
Research Area	Others
Target Names	Major allergen Bla g 1.0101
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.