



Allergen Bla g 4 Antibody

Product Code	CSB-PA19857A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54962
Immunogen	Recombinant Blattella germanica Allergen Bla g 4 protein (13-182AA)
Raised In	Rabbit
Species Reactivity	Blattella germanica
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Allergen Bla g 4 (Allergen Bla g IV) (allergen Bla g 4) (Fragment)
Immunogen Species	Blattella germanica (German cockroach) (Blatta germanica)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.