




Major pollen allergen Aln g 1 Antibody

Product Code	CSB-PA19897A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38948
Immunogen	Recombinant Alnus glutinosa Major pollen allergen Aln g 1 protein (2-160AA)
Raised In	Rabbit
Species Reactivity	Alnus glutinosa
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Aln g 1 (Allergen Aln g I) (allergen Aln g 1)
Immunogen Species	Alnus glutinosa
Research Area	Others

Image



Western Blot

Positive WB detected in Recombinant protein
 All lanes: Major pollen allergen Aln g 1 antibody at 3µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 predicted band size: 34 kDa
 observed band size: 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.