



Major pollen allergen Aln g 1 Antibody, HRP conjugated

Product Code	CSB-PA19897B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38948
Immunogen	Recombinant Alnus glutinosa Major pollen allergen Aln g 1 protein (2-160AA)
Raised In	Rabbit
Species Reactivity	Alnus glutinosa
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Aln g 1 (Allergen Aln g I) (allergen Aln g 1)
Immunogen Species	Alnus glutinosa (European alder) (Betula alnus var. glutinosa)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.