




BETV4 Antibody

Product Code	CSB-PA19907A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39419
Immunogen	Recombinant Betula pendula Polcalcin Bet v 4 protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Betula pendula
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polcalcin Bet v 4 (Calcium-binding pollen allergen Bet v 4) (allergen Bet v 4), BETV4
Immunogen Species	Betula pendula
Research Area	Others
Target Names	BETV4

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: BETV4 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 26 kDa

observed band size: 26 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.