



BETV4 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA19907B0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q39419 |
| Immunogen | Recombinant Betula pendula Polcalcin Bet v 4 protein (1-85AA) |
| Raised In | Rabbit |
| Species Reactivity | Betula pendula |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Polcalcin Bet v 4 (Calcium-binding pollen allergen Bet v 4) (allergen Bet v 4), BETV4 |
| Immunogen Species | Betula pendula (European white birch) (Betula verrucosa) |
| Research Area | Others |
| Target Names | BETV4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |