



BETV4 Antibody, Biotin conjugated

Product Code	CSB-PA19907D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39419
Immunogen	Recombinant Betula pendula Polcalcin Bet v 4 protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Betula pendula
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polcalcin Bet v 4 (Calcium-binding pollen allergen Bet v 4) (allergen Bet v 4), BETV4
Immunogen Species	Betula pendula (European white birch) (Betula verrucosa)
Research Area	Others
Target Names	BETV4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.