



OLE9 Antibody, HRP conjugated

Product Code	CSB-PA19957B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q94G86
Immunogen	Recombinant Olea europaea Glucan endo-1,3-beta-D-glucosidase protein (360-460AA)
Raised In	Rabbit
Species Reactivity	Olea europaea
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucan endo-1,3-beta-D-glucosidase (EC 3.2.1.39) (Major pollen allergen Ole e 9) (allergen Ole e 9), OLE9
Immunogen Species	Olea europaea (Common olive)
Research Area	Others
Target Names	OLE9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.