



egfra Antibody, FITC conjugated

Product Code	CSB-PA20899C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F1RA48
Immunogen	Recombinant Zebrafish Epidermal growth factor receptor protein (21-389AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	egfra
Immunogen Species	Zebrafish
Research Area	Others
Target Names	egfra
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.