





ALK Antibody

Product Code	CSB-PA21109A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM73
Immunogen	Recombinant Human ALK tyrosine kinase receptor protein (1402-1604AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ALK tyrosine kinase receptor (EC 2.7.10.1) (Anaplastic lymphoma kinase) (CD antigen CD246), ALK
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ALK

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA21109A0Rb at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA21109A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.