





PPP2R3A Antibody

Product Code	CSB-PA21129A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06190
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit alpha protein (256-508AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit alpha (PP2A subunit B isoform PR72/PR130) (PP2A subunit B isoform R3 isoform) (PP2A subunit B isoforms B"-PR72/PR130) (PP2A subunit B isoforms B72/B130) (Serine/threonine-protein phosphatase 2A 72/130 kDa regulatory subunit B), PPP2R3A, PPP2R3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPP2R3A

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA21129A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA21129A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.