



PPP2R3A Antibody, Biotin conjugated

Product Code	CSB-PA21129D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06190
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit alpha protein (256-508AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit alpha (PP2A subunit B isoform PR72/PR130) (PP2A subunit B isoform R3 isoform) (PP2A subunit B isoforms B"-PR72/PR130) (PP2A subunit B isoforms B72/B130) (Serine/threonine-protein phosphatase 2A 72/130 kDa regulatory subunit B), PPP2R3A, PPP2R3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPP2R3A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.