



VII Antibody

Product Code	CSB-PA302477LA01ECY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69535
Immunogen	Recombinant Enterobacteria phage M13 Tail virion protein G7P protein (1-33AA)
Raised In	Rabbit
Species Reactivity	Enterobacteria phage M13
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tail virion protein G7P (Coat protein C, polypeptide I) (Gene 7 protein) (G7P), VII
Immunogen Species	Enterobacteria phage M13 (Bacteriophage M13)
Research Area	Others
Target Names	VII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.