



## VII Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA302477LD01ECY
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P69535
<b>Immunogen</b>	Recombinant Enterobacteria phage M13 Tail virion protein G7P protein (1-33AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Enterobacteria phage M13
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Tail virion protein G7P (Coat protein C, polypeptide I) (Gene 7 protein) (G7P), VII
<b>Immunogen Species</b>	Enterobacteria phage M13 (Bacteriophage M13)
<b>Research Area</b>	Others
<b>Target Names</b>	VII