




GP Antibody

Product Code	CSB-PA308957LA01ZAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P87670
Immunogen	Recombinant Zaire ebolavirus Pre-small/secreted glycoprotein protein (33-364AA)
Raised In	Rabbit
Species Reactivity	Zaire ebolavirus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-small/secreted glycoprotein (pre-sGP) [Cleaved into: Small/secreted glycoprotein (sGP); Delta-peptide], GP
Immunogen Species	Zaire ebolavirus
Research Area	Others
Target Names	GP

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: GP antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

predicted band size: 42 kDa

observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.