



# GP Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA308957LD01ZAA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P87670
<b>Immunogen</b>	Recombinant Zaire ebolavirus Pre-small/secreted glycoprotein protein (33-364AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Zaire ebolavirus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pre-small/secreted glycoprotein (pre-sGP) [Cleaved into: Small/secreted glycoprotein (sGP); Delta-peptide], GP
<b>Immunogen Species</b>	Zaire ebolavirus
<b>Research Area</b>	Others
<b>Target Names</b>	GP