



rgpB Antibody, HRP conjugated

Product Code	CSB-PA310587LB01EYA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P95493
Immunogen	Recombinant Porphyromonas gingivalis Gingipain R2 protein (230-473AA)
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gingipain R2 (EC 3.4.22.37) (Arg-gingipain) (Gingipain 2) (RGP-2), rgpB, prtRII rgp2
Immunogen Species	Porphyromonas gingivalis
Research Area	Others
Target Names	rgpB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.