



CDH-1 Antibody, HRP conjugated

Product Code	CSB-PA312474YB01EUK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01738
Immunogen	Recombinant Phanerochaete chrysosporium Cellobiose dehydrogenase protein (19-208AA)
Raised In	Rabbit
Species Reactivity	Phanerochaete chrysosporium
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellobiose dehydrogenase (CDH) (EC 1.1.99.18) (Cellobiose-quinone oxidoreductase), CDH-1; CDH-2
Immunogen Species	Phanerochaete chrysosporium (White-rot fungus) (Sporotrichum pruinosum)
Research Area	Others
Target Names	CDH-1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.