



pfo Antibody, HRP conjugated

Product Code	CSB-PA314820HB01CMB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C2E9
Immunogen	Recombinant Clostridium perfringens Perfringolysin O protein (29-500AA)
Raised In	Rabbit
Species Reactivity	Clostridium perfringens
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Perfringolysin O (Theta-toxin) (Thiol-activated cytolysin), pfo, pfoA pfoR
Immunogen Species	Clostridium perfringens (strain 13 / Type A)
Research Area	Others
Target Names	pfo
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.