



Defensin-like protein 1 Antibody, HRP conjugated

Product Code	CSB-PA316614LB01DAE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C8Y4
Immunogen	Recombinant Dahlia merckii Defensin-like protein 1 protein (1-50AA)
Raised In	Rabbit
Species Reactivity	Dahlia merckii
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Defensin-like protein 1 (Cysteine-rich antimicrobial protein 1) (Defensin AMP1) (DmAMP1)
Immunogen Species	Dahlia merckii (Bedding dahlia)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.