



N Antibody, FITC conjugated

Product Code	CSB-PA318492LC01JAM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14239
Immunogen	Recombinant Junin mammarenavirus Nucleoprotein protein (1-564AA)
Raised In	Rabbit
Species Reactivity	Junin mammarenavirus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), N
Immunogen Species	Junin mammarenavirus
Research Area	Others
Target Names	N
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.