



# mrkA Antibody, Biotin conjugated

Product Code	CSB-PA318939LD01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12267
Immunogen	Recombinant Klebsiella pneumoniae Fimbrial subunit type 3 protein (23-202AA)
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fimbrial subunit type 3, mrkA
Immunogen Species	Klebsiella pneumoniae
Research Area	Others
Target Names	mrkA