



mrkA Antibody, Biotin conjugated

Product Code	CSB-PA318939LD01KBG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12267
Immunogen	Recombinant Klebsiella pneumoniae Fimbrial subunit type 3 protein (23-202AA)
Raised In	Rabbit
Species Reactivity	Klebsiella pneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fimbrial subunit type 3, mrkA
Immunogen Species	Klebsiella pneumoniae
Research Area	Others
Target Names	mrkA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.