



Sporamin A Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA320416HD01INB |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P14715 |
| Immunogen | Recombinant Ipomoea batatas Sporamin A protein (24-219AA) |
| Raised In | Rabbit |
| Species Reactivity | Ipomoea batatas |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Sporamin A (Clone PIM0335) |
| Immunogen Species | Ipomoea batatas (Sweet potato) (Convolvulus batatas) |
| Research Area | Others |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |