



bla Antibody

Product Code	CSB-PA321251LA01EZW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14489
Immunogen	Recombinant Pseudomonas aeruginosa Beta-lactamase OXA-10 protein (21-266AA)
Raised In	Rabbit
Species Reactivity	Pseudomonas aeruginosa
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-lactamase OXA-10 (EC 3.5.2.6) (Beta-lactamase PSE-2), bla, oxa10 pse2
Immunogen Species	Pseudomonas aeruginosa
Research Area	Others
Target Names	bla
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.