



BETVIA Antibody, HRP conjugated

Product Code	CSB-PA322213FHB01BSS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15494
Immunogen	Recombinant Betula pendula Major pollen allergen Bet v 1-A protein (2-160AA)
Raised In	Rabbit
Species Reactivity	Betula pendula
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major pollen allergen Bet v 1-A (Allergen Bet v I-A) (allergen Bet v 1-A), BETVIA, BETVI
Immunogen Species	Betula pendula (European white birch) (Betula verrucosa)
Research Area	Others
Target Names	BETVIA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.