




D8L Antibody

Product Code	CSB-PA322653LA01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20508
Immunogen	Recombinant Vaccinia virus Cell surface-binding protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface-binding protein (Carbonic anhydrase homolog)
Immunogen Species	Vaccinia virus
Research Area	Others
Target Names	D8L

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: D8L antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa