



D8L Antibody, FITC conjugated

Product Code	CSB-PA322653LC01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20508
Immunogen	Recombinant Vaccinia virus Cell surface-binding protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface-binding protein (Carbonic anhydrase homolog)
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Research Area	Others
Target Names	D8L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.