



D8L Antibody, Biotin conjugated

Product Code	CSB-PA322653LD01VAA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20508
Immunogen	Recombinant Vaccinia virus Cell surface-binding protein (1-304AA)
Raised In	Rabbit
Species Reactivity	Vaccinia virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell surface-binding protein (Carbonic anhydrase homolog)
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Research Area	Others
Target Names	D8L