



p37 Antibody, HRP conjugated

Product Code	CSB-PA323355HB01MLR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15363
Immunogen	Recombinant Mycoplasma hyorhinis High affinity transport system protein p37 protein (24-403AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma hyorhinis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity transport system protein p37, p37
Immunogen Species	Mycoplasma hyorhinis
Research Area	Others
Target Names	p37