



p37 Antibody, Biotin conjugated

Product Code	CSB-PA323355HD01MLR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15363
Immunogen	Recombinant Mycoplasma hyorhinis High affinity transport system protein p37 protein (24-403AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma hyorhinis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity transport system protein p37, p37
Immunogen Species	Mycoplasma hyorhinis
Research Area	Others
Target Names	p37
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.