



p37 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA323355HD01MLR |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P15363 |
| Immunogen | Recombinant Mycoplasma hyorhinis High affinity transport system protein p37 protein (24-403AA) |
| Raised In | Rabbit |
| Species Reactivity | Mycoplasma hyorhinis |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | High affinity transport system protein p37, p37 |
| Immunogen Species | Mycoplasma hyorhinis |
| Research Area | Others |
| Target Names | p37 |